

Quick Start Guide

This is a shorter version of the user guide that helps the users to get familiar with the platform quickly.

- Steps to perform your first analysis
- Customisation

Steps to perform your first analysis

Step 1 : Upload samples and other required files through "Uploader". Each file will go through the validation process. Check the [DataStore](#) to verify the status of the uploaded files.

Step 2 : Create your first project using [Projects](#) from the left menu.

Step 3 : Inside the project, create a [Sample Set](#) from the [DataStore](#).

Step 4: Create a new pipeline or copy a preconfigured pipeline from the PipelineLibrary ([My Pipelines](#)).





Step 5: Execute the pipeline by clicking the [Initialize Pipeline](#) button. Fill in the required fields and click "Execute Pipeline".

Step 6: Monitor pipeline execution progress through [Jobs](#).

Step 7: After the run completion, [Report](#) is ready to reviewed.

Customisation

Adding a New Organism

Click on settings icon  from the right upper corner. Select  Organism from the left menu bar. Create a new organism by clicking  icon. Fill in the required information and click  icon to save it.