

# Quick Start Guide

This is a shorter version of the user guide that helps the users to get familiar with the platform quickly.

- [Steps to perform your first analysis](#)
- [Customisation](#)

# Steps to perform your first analysis

**Step 1 :** Upload samples and other required files through "[Uploader](#)". Each file will go through the validation process. Check the [DataStore](#) to verify the status of the uploaded files.

**Step 2 :** Create your first project using [Projects](#) from the left menu.

**Step 3 :** Inside the project, create a [Sample Set](#) from the [DataStore](#).

**Step 4:** Create a new pipeline or copy a preconfigured pipeline from the PipelineLibrary ([My Pipelines](#)).





**Step 5:** Execute the pipeline by clicking the [Initialize Pipeline](#) button. Fill in the required fields and click "Execute Pipeline".

**Step 6:** Monitor pipeline execution progress through [Jobs](#).

**Step 7:** After the run completion, [Report](#) is ready to reviewed.

# Customisation

## Adding a New Organism

Click on settings icon  from the right upper corner. Select  Organism from the left menu bar. Create a new organism by clicking  icon. Fill in the required information and click  icon to save it.